EN

DA

DE

EL

ES

FI

FR

IT

JA

NL

NO

PL

PT

SV

TR

## INSTRUCTIONS FOR USE

# ISODRY®

Dental Isolation System

TABLE OF CONTENTS  
INDHOLDSFORTEGNELSE  
INHALTSVERZEICHNIS  
ΠΙΝΑΚΑΣ ΠΕΡΙΕΧΟΜΕΝΩΝ  
CONTENIDO  
SISÄLLYSLUETTELO  
TABLE DES MATIÈRES  
INDICE  
目次  
INHOUDSOPGAVE  
INNHOLDSFORTEGNELSE  
SPIS TREŚCI  
ÍNDICE  
INNEHÅLLSFÖRTECKNING  
İNDEKİLER

**ISODRY®**

Isodry® Dental Isolation System: Instructions for Use

EN

Isodry®-tandisolationssystem: Brugsanvisning

DA

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[REF](#) EID0103 (EU)

[REF](#) EID0104 (AUS)

This system is intended for use with Mouthpieces developed by Zyris, Inc. These Mouthpieces are available in a wide-range of sizes for effective isolation.

Each Isodry Dental Isolation System is comprised of three main components: 1) Titanium Control Base, 2) Vacuum Pipe, 3) Vacuum Hose.

#### Titanium Control Base

[REF](#) PID0101

The Titanium Control Base includes two Vacuum Levers (P/N: PIL0205) to adjust suction in the upper and lower areas of the Mouthpiece.

#### Vacuum Pipe

[REF](#) PID0103

The Vacuum Pipe channels suction to the Mouthpiece. The Mouthpiece attaches onto the Vacuum Pipe.

#### Vacuum Hose

[REF](#) PID0102 (USA)

(1.82 meters / 72 inches)

[REF](#) PID0104 (EU)

(1.49 meters / 59 inches)

[REF](#) PID0102 (NZ/AUS)

(1.82 meters / 72 inches)

The Vacuum Hose delivers vacuum suction from the

vacuum canister to the Control Head. The Vacuum Hose connects to the vacuum canister.

#### Suction Line Plug

[REF](#) PIL0206

#### O-ring Maintenance Kit

[REF](#) AIL0206 (USA/NZ/AUS)

[REF](#) AIL0223 (EU)

#### 3/8" Y-adapter Kit

[REF](#) AIL0207 (USA/AUS)

#### 9.5 mm Y-adapter Kit

[REF](#) AIL0222 (EU/AUS)

#### Double End Tube Brush

[REF](#) AIL0204

#### Mouthpiece Lubricant

[REF](#) AIL0202

#### Barrier Sleeves

[REF](#) AIL0203 (non-sterile)

Sold Separately

Non-sterile Control Head Barrier Sleeves are available for additional protection from cross-contamination. Use of Barrier Sleeves does not eliminate the need to follow disinfecting and sterilization recommendations for the Control Head and Vacuum Pipe.

### Mouthpieces

REF CILXXOX (non-sterile)

Sold Separately

Mouthpieces are single use only and available in various sizes. The Mouthpieces are constructed of a soft, flexible polymer and are not made with natural rubber latex. They do not contain phthalates or BPAs. For more information on sizes, go to Zyris.com.

### Intended Use

Isodry is a dental isolation system for clinical use by dental professionals. The system is designed to isolate two quadrants at once, simultaneously providing tongue and cheek retraction and continuous suction. The retraction and continuous suction are the important functions that control the oral environment, reducing contamination and improving bond strength while protecting the airway and soft tissue.

### Contraindications

No known contraindications. Do not use without reviewing how to select the correct Mouthpiece size for each patient. Do not use if unable to position the Mouthpiece correctly.

### Side Effects

The use of the Mouthpiece may trigger a pharyngeal reflex (laryngeal spasm) in hypersensitive patients, or if it is not placed correctly within the oral cavity.

### How Supplied

#### Isodry

The Isodry Dental Isolation System contains: Titanium Control Base, 6 Vacuum Pipes, Vacuum Hose, and Instructions for Use.

Accessories: Suction Line Plug, Mouthpiece Lubricant, Y-adaptor Kit, Double End Tube Brush, and O-ring Maintenance Kit.

### Mouthpieces

Mouthpieces (sold separately) are supplied in convenient operatory boxes. Each Mouthpiece is individually packaged.

## Installation

Isodry is designed to install easily in most dental operatories (Figure 1). Installation and usage instructions, including a video demonstration, are available at [Zyris.com](http://Zyris.com).

**Step 1** Connect the Vacuum Hose to the vacuum canister by fitting the end of the Vacuum Hose onto an auxiliary HVE port on the vacuum canister. If an auxiliary HVE port is not available, use the provided Y-adapter, so that the standard HVE hose may remain connected.


Figure 1. Installation

*NOTE: For countries outside the USA, installation may need to be completed by a qualified service technician. Please consult your in-country representative for additional information.*

**Step 2** Connect the Control Head to the head of the Vacuum Hose. Once connected, test the dual Vacuum Levers to confirm suction functionality.

**Step 3** Place the Isodry in an open slot on the toolbar when not in use.

## How to Operate

### Using the Dual Vacuum Levers

Dual Vacuum Levers adjust suction intensity in both the upper and lower quadrants. To adjust, rotate the Left and Right Vacuum Levers to the desired suction intensity for each quadrant. To shut off suction, rotate the Levers to the most downward position.

### Selecting the Mouthpiece

To select the correct Mouthpiece size, measure the patient's interincisal opening. If between sizes, you may select the smaller size.

Mouthpiece Size	Approx. Interincisal Opening
Pediatric (PED)	< 30 mm
Extra Small (XSM)	< 30 mm
Small (SML)	30 mm
Medium (MED)	30-45 mm
Medium Deep Vestibule (MDV)	30-45 mm
Large (LRG)	> 45 mm

Table 1. Mouthpiece Sizes

### The Finger Method

The finger method may also be used to select the appropriate Mouthpiece size (Figure 2). For additional assistance, watch a video demonstration online at Zyris.com.


Figure 2. The Finger Method


Figure 3. Mouthpiece

### Attaching the Mouthpiece

Slide the Mouthpiece onto the Vacuum Pipe until the two locking bars engage. Water or mouthpiece lubricant may be used to lubricate the inside surface of the Mouthpiece sleeve to ease attachment.

### Placing the Mouthpiece

Prior to placing the Mouthpiece within the patient's oral cavity, moisten the Mouthpiece and the patient's lips.

#### Step 1

Hold the sides of the Control Head with your index finger and thumb. With your other hand, fold the Cheek Shield onto the Tongue Retractor. Gently slide the folded

mouthpiece into the buccal vestibule on the side to receive treatment. Angle the mouthpiece to allow the lower edge of the Tongue Retractor to move along the buccal edge of the teeth.

### Step 2

Move the bite block onto the occlusal surface of the teeth, just distal to the mandibular cuspid. Instruct the patient to "rest gently" on the bite block to secure.


Figure 4, 5, & 6. Mouthpiece Placement

### Step 3

Place the Cheek Shield into the buccal vestibule. Move the Isthmus onto the retromolar pad behind the maxillary tuberosity. Adjust the Tongue Retractor in the lingual vestibule as needed. Move the bite block distally to provide more vertical working room.

## Disassembly of the Control Head

### Step 1

Press the release button (Figure 7) on the underside of the Titanium Control Base to release and disconnect the Vacuum Pipe.

### Step 2

Holding the head of the Vacuum Hose, pull to remove the Titanium Control Base. Do not twist.


Figure 7. Control Head

## Disinfection, Cleaning, and Sterilization

### Surface Disinfection of the Titanium Control Base

Zyris recommends performing a surface disinfection of the Titanium Control Base between each patient. To disinfect the Titanium Control Base, first remove the Vacuum Pipe.

If the Titanium Control Base is visibly soiled with biomatter or blood, use an intermediate level disinfectant wipe with a claim to inactivate *Mycobacterium tuberculosis*. Otherwise, use a hospital grade disinfectant and follow the manufacturer's instructions for contact time.

### Surface Disinfection of the Vacuum Hose

Zyris recommends performing a surface disinfection of the Vacuum Hose between each patient.

Use a hospital grade disinfectant per manufacturer's instructions. If the components are visibly soiled with biomatter or blood, use an intermediate level disinfectant wipe with a claim to inactivate *Mycobacterium tuberculosis*. Barrier Sleeves may also be used and changed between patients.

### Cleaning and Sterilization of the Vacuum Pipe

Zyris recommends sterilization of the Vacuum Pipe between each patient.

Remove biomatter and other debris from the suction channels using the provided Double End Tube Brush with a neutral pH (6.5-7.5) enzymatic cleaning solution. Do not use solution to hold parts. After cleaning, rinse with tap water, dry, and place components in an autoclave pouch to be sterilized in a steam autoclave. Use validated sterilization cycles (Table 2).

### Cleaning and Sterilization of the Titanium Control Base

Zyris recommends following normal infection control protocol for dental devices. Users may sterilize the Titanium Control Base at their discretion.

Remove Vacuum Levers from the Titanium Control Base. Remove biomatter and other debris from the suction channels using the provided Double End Tube Brush with a neutral pH (6.5-7.5) enzymatic cleaning solution. Do not hold parts in cleaning solution or use ultrasonic cleaners. After cleaning, rinse with tap water, dry, and place components in an autoclave pouch to be sterilized in a steam autoclave. Use validated sterilization cycles (Table 2).

After sterilization, lubricate the Vacuum Lever O-rings with the Dow 111 O-ring lubricant

provided in the O-ring Maintenance kit.

#### Validated Sterilization Cycles

Cycle	Temperature	Min. Exposure Time	Min. Drying Time
Standard	121 °C / 250 °F	30 min	30 min
Standard	132 °C / 270 °F	15 min	15 min
Prevacuum	132 °C / 270 °F	4 min	20 min
Prevacuum	134 °C / 273 °F	3 min	20 min

Table 2. Validated Sterilization Cycles

## Maintenance

### Titanium Control Base

The Vacuum Lever O-rings require weekly lubrication to operate at peak efficiency. Remove the Titanium Control Base from the Vacuum Hose. To remove the Vacuum Levers, push the levers to the "off" position and keep pushing gently until each lever pops out, or gently pry them out. Coat the O-ring by applying a small amount of the Dow 111 O-ring Lubricant provided in the O-ring Maintenance Kit. Apply an additional amount to the Titanium Control Base retainer rings as well. Re-attach the Vacuum Levers to the Titanium Control Base (Figure 8 and 9).

Figure 8. Control Base

Figure 9. Vacuum Levers

### Vacuum Hose

Clean the inside of the Vacuum Hose regularly according to frequency of use to prevent the accumulation of biofilm. Start by vacuuming water through the Vacuum Hose. Then vacuum a dental hose cleaner through the hose per manufacturer's instructions. The vacuum hose cleaner should be non-foaming, deodorizing, anti-microbial, and able to dissolve organic debris. Most enzymatic cleaners with pH between 6 and 8 are acceptable. Do not use disinfectants as a vacuum hose cleaner. If further disinfection is necessary, consult your vacuum manufacturer for recommendations.

To lubricate the Vacuum Hose O-ring, remove the Titanium Control Base from the Vacuum Hose. The black O-ring is located approximately 13 mm (½ inch) down the inside of the blue fitting at the head of the

Vacuum Hose. (See Figure 10) Swab the O-ring with the Dow 111 O-ring lubricant provided in the O-ring Maintenance Kit.

Figure 10. Vacuum Hose

## Warning

**WARNING:** Mouthpieces are single use only. Reuse can result in the transmission of infectious pathogens. Mouthpieces must be discarded after use. Mouthpieces cannot be safely resterilized after use.

**WARNING:** Do not use if the individual Mouthpiece package has been previously damaged or opened. Non-sterile Mouthpieces are individually packaged in polyethylene bags.

**WARNING:** Be sure to use established protocols and patient monitoring when using the Mouthpiece in conjunction with sedation dentistry. Patient management by a certified anesthesiologist is strongly recommended.

## Precautions

**CAUTION:** When performing dental procedures requiring high levels of blood evacuation, it is recommended that the Mouthpiece be intermittently flushed with water to prevent coagulation in the Mouthpiece to avoid loss of suction.

**CAUTION:** Do not autoclave Titanium Control Base at temperatures above 135 °C; higher temperatures could damage components and void your product warranty.

**CAUTION:** Do not hold Titanium Control Base in cleaning solution or use ultrasonic cleaners.

**CAUTION:** Failure to maintain properly lubricated O-rings can result in damage to the Isodry that is not covered by product warranty.

**CAUTION:** With regular use, cleaning, and maintenance, the device is designed to last for years. Check the device for wear prior to use to make sure that it is adequate for use.

## Guide to Symbols

Catalog number	Non-sterile	Caution
Date of manufacture	Do not re-use	Do not use if package is damaged
Consult Instructions for Use	Manufacturer	Storage Conditions

## Limited 1-Year Warranty

Zyris warrants this product to be free from defects in material and workmanship for a period of one year from the date of original purchase. If the product exhibits such a defect, Zyris will, at its option, repair or replace it without costs for parts and labor. The product must be returned by the customer to Zyris in accordance with current Return Material Authorization procedures. Proof of date of original purchase may be required. This warranty does not cover finishes or normal wear, nor does it cover damage resulting from accident, misuse, dirt, tampering, unreasonable use, failure to provide reasonable and necessary maintenance, service performed or attempted by unauthorized individuals or service agencies, or devices that have been otherwise modified. All implied warranties, including any implied warranty or merchantability or fitness for any particular purpose, are limited in duration to one year from the date of receipt of the product. In no event will Zyris be responsible for consequential damages resulting from the use of this product.

Instructions for Use may also be found at: [Zyris.com](http://Zyris.com).

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All rights reserved.
© 2018 Zyris, Inc.