

HOOPLA™

Sizing Guide

The ultimate dental retractor that's so simple and comfortable to use your patient could insert it themselves in seconds!

Two Sizes Available

Regular: Fits most patients

[CIR0204] 40 Pack Regular
[CIR0208] 80 Pack Regular

Small: Fits small children and adults with a small mouth

[CIR0304] 40 Pack Small
[CIR0308] 80 Pack Small

Dimensions

	Regular	Small
Length	154 mm	120 mm
Width	75 mm	58 mm
Lip Tray Width	54 mm	42 mm
Folded Length	73 mm	56 mm
Intraoral Depth	34 mm	27 mm

